



PROCEDURE PACKS CLINICAL CONSUMABLES SINGLE-USE INSTRUMENTS

www.mddiglobal.com

mddi | 
STERILE | SINGLE USE

PROCEDURE PACKS



Utilise our
drape kits
during the
following
procedures

Implant
Surgery

Minor Oral
Surgery

Surgical
Extractions

Periodontal
Surgery

Our kits contain high-quality components, thoughtfully constructed through market research and direct collaboration with dental practitioners. Featuring a selection of drapes and gowns of varying sizes and absorbent qualities as well as an array of clinical consumables such as facemasks, suction tubing and x-ray gauze.

PP070ST-1 Standard Implant and Oral Surgery Procedure Pack

3 packs per box | 2 boxes per carton



*Please note images are not to scale

CONTENTS

2 x PA001NS

Surgical Gown 60gsm Sterile
AAMI Level 2

1 x SD015NS

Surgical Patient Drape Impervious
U-shape (180cm x 270cm) +
Adhesive + Tubeholders

1 x SD019NS

Mayo Stand Cover (140cm x 60cm)

1 x SD020NS

Head Bar Drape with Adhesive
(80cm x 100cm)

1 x SD021NS

Trolley/Back Table Cover
(180cm x 140cm)

4 x CC009NS

Hand Towel (40cm x 40cm)

3 x CC022NS

Adhesive Film Blue (20cm x 20cm)

1 x CC023NS

Fine Aspiration Handle and Tip
with Control Vent (2.8mmID)

1 x CC024NS

Blue Adaptor

10 x CC031NS

Gauze X-Ray Detectable
8ply (10cm x 10cm)

1 x CC033NS

Saliva Ejector

1 x CC034NS

Tubing Adaptor/Connector

1 x CC035NS

Suction Tube with Connectors (3m)

2 x CC075NS

Quick Drill Sleeve with Applicator
(150cm x 9.5cm)



PP071ST-1 Economy Implant and Oral Surgery Procedure Pack

8 packs per box | 2 boxes per carton



SD013NS



CC075NS

SD011NS




CC036NS

CC022NS



CC031NS

SD022NS



SD016NS



CONTENTS

*Please note images are not to scale

1 x SD011NS

Surgical Drape Impervious (90cm x 120cm)

1 x SD013NS

Surgical Drape Impervious (75cm x 90cm)

1 x SD016NS

3/4 Patient Drape/U Drape (100cm x 150cm)
+ Adhesive + Tube Holders

1 x SD022NS

Trolley/Back Table Cover (140cm x 140cm)

4 x CC022NS

Adhesive Film Blue (20cm x 20cm)

16 x CC031NS

Gauze X-Ray Detectable 8ply (10cm x 10cm)

1 x CC036NS

Biohazard Waste Bag Yellow (42cm x 45cm)

2 x CC075NS

Quick Drill Sleeve with Applicator (150cm x 9.5cm)

PP074ST-1 Complete Implant and Oral Surgery Procedure Pack

3 packs per box | 2 boxes per carton



CONTENTS

*Please note images are not to scale

1 x CC024NS

Blue Adaptor

1 x SD011NS

Surgical Drape Impervious
(90cm x 120cm)

1 x SD013NS

Surgical Drape Impervious
(75cm x 90cm)

1 x SD021NS

Trolley/Back Table Cover
(180cm x 140cm)

3 x CC022NS

Adhesive Film Blue
(20cm x 20cm)

3 x CC012NS

Adhesive Film Clear
(5cm x 80cm)

1 x CC023NS

Fine Aspiration Handle and Tip
with Control Vent (2.8mmID)

3 x CC027NS

Drill Sleeve with Ties Clear
(120cm)

1 x CC034NS

Tubing Adaptor/Connector

1 x CC035NS

Suction Tube with Connectors
(3m)

1 x CC036NS

Biohazard Waste Bag Yellow
(42cm x 45cm)

1 x CC037NS

Fine Tip Yankauer Suction
Handle with Vent (3.2mmID)

1 x CC094NS

Blue Plastic Tweezers

1 x CC028NS

Bouffant Style
Equipment Cover Clear
(45cm)

1 x CC029NS

Gallipot Clear (60ml)

1 x CC030NS

Kidney Dish Clear
(700ml)

16 x CC031NS

Gauze X-Ray Detectable 8ply
(10cm x 10cm)

4 x CC009NS

Hand Towel (40cm x 40cm)

1 x SD016NS

3/4 Patient Drape/U Drape
(100cm x 150cm) with
Adhesive + Tube Holders

2 x PA001NS

Surgical Gown 60gsm
Sterile AAMI Level 2

2 x PA0014NS

Type IIR Face Mask
with Ear Loops

1 x PA004NS

Type IIR Face Mask with Ties

2 x PA005NS

Bouffant Cap/Head
Cover with Elastic

1 x PA006NS

Tie Back Head Cover





CLINICAL CONSUMABLES

In response to critical demand MDDI has developed a cost effective and innovative range of single-use clinical consumables suitable for use within both dental and medical environments.

SURGICAL DRAPES



Surgical Gown Bundle

PA073ST-1

30 gown bundles per carton

1 x Surgeon Gown (125cm x 155cm)

60gsm AAMI Level 2

2 x Hand Towel (40cm x 40cm)



Sterile Surgical Drape Impervious (90cm x 120cm)

SD011ST-1

18 drapes per box

8 boxes per carton



Sterile Surgical Drape Impervious (75cm x 90cm)

SD013ST-1

25 drapes per box

8 boxes per carton




Sterile 3/4 Patient Drape/U Drape with Adhesive + Tube Holders (100cm x 150cm)

SD016ST-1

10 drapes per box

8 boxes per carton





Aspirator Kit

CC072ST-1

8 kits per box | 8 boxes per carton

1 x Suction Tube with Connectors (3m)


1 x Fine Aspiration Handle and Tip with Control Vent (2.8mmID)

Irrigation Tubing - Set A

CC101ST-1

20 sets per box | 8 boxes per carton

Compatible with NSK – Surgic XT/Surgic XT PLUS/Variosurg, Biomet 3I – Osseocision, KYOCERA – Implantor Neo. ASEPTICO – AEU 7000/AEU 6000 motors.



Irrigation Tubing - Set B

CC102ST-1

20 sets per box | 8 boxes per carton


Compatible with Nobel Biocare80, W&H Implantmed, Implant Med new design 2005, Implantmed SI 923, Astratech Implant SI 923, Straumann, Osseoset 100, Osseoset 200, Osseocare.



Surgical Face Mask with Elasticated Ear Loops - Level 3

PA014NS-1

50 masks per box
40 boxes per carton



Our face masks are CE certified conforming to AS 4381:2015

- Non-sterile • Latex free • Earloops • 3-ply non-woven fabric with 25gsm filtration
- Fluid resistance 160mmHg • Delta-P <5.0 • PFE >99% @ 0.1µm • BFE >98%

BARRIERS



Drill Sleeve Elasticised Clear on Card (120cm x 7.5cm) Pkt1

CC126ST-1

50 individual sleeves packaged as 50 packs of 1
8 boxes per carton



Drill Sleeve Elasticised Clear on Card (120cm) Pkt2

CC026ST-1

60 individual sleeves packaged as 30 packs of 2
8 boxes per carton



Quick Drill Sleeve with Applicator (150cm x 9.5cm) Pkt1

CC175ST-1

24 individual sleeves packaged as 24 packs of 1
8 boxes per carton



Quick Drill Sleeve with Applicator Ring (150cm x 9.5cm) Pkt2

CC075ST-1

24 individual sleeves packaged as 12 packs of 2
8 boxes per carton



Adhesive Film Blue (20cm x 20cm) Pkt 3

CC022ST-1


50 packs per box | 6 boxes per carton



IMPRESSION TRAYS

Single-use plastic impression trays with integrated handle. Constructed from rigid, recyclable ABS for accurate, comfortable and distortion free impressions.

Non-Sterile. 25 trays per box.



Impression Tray Small Upper
IT02SMU-1



Impression Tray Small Lower
IT01SML-1



Impression Tray Medium Upper
IT04MBU-1



Impression Tray Medium Lower
IT03MBL-1



Impression Tray Large Upper
IT06LGU-1



Impression Tray Large Lower
IT05LGL-1



Impression Tray Extra Large Upper
IT08XLU-1



Impression Tray Extra Large Lower
IT07XLL-1





SINGLE-USE INSTRUMENTS

Our Instruments are:

- Quick and easy to set up and use
- Prewrapped in a range of convenient kits and instruments
- Sterile at the point of use
- Convenient to dispose of in sharps container
- Reducing the risk of sharps injuries
- Sharp where needed, ensuring the highest quality of treatment for every procedure
- Implemented to significantly reduce the risk of cross-contamination
- GS1 compliant with removable tracking labels
- CE certified and TGA registered
- 5 year expiry from point of manufacture

MDDI Can Increase Revenue by:

- Minimise decontamination cycles per shift
 - reduces expenditure on consumables, test kits and utilities
- Putting the RDN back chairside to assist the dentist and patient
- Increasing the number of patients seen in each clinical session
- Improving clinical efficiencies allowing for longer treatment plans and procedures
- Improving the patient journey



STERILE, SINGLE-USE DIAGNOSTIC KITS



Our MDDI diagnostic kits were adopted into the following sectors more than 15 years ago. We have continued to expand our range due to the ongoing popularity of our instruments.

School
Clinics

Government
Hospitals

Flying Doctors
Service

Mobile
Dental Units

Our *diagnostic instruments* are high quality, precision made instruments. The probes and London College Tweezers are manufactured from stainless steel and all of our handles are recyclable ABS plastic.

INDIVIDUAL DIAGNOSTIC

Front Face Mirror with Percussion Tip

SI053ST-1 - 180 per box




London College Tweezers

SI182ST-1 - 200 per box



Tweezers

SI081ST-1 - 200 per box





Periodontal 2 Piece Examination Kit with CPITN-C Probe

SI952ST-1 - 100 kits per box

1 x Front Face Mirror with Percussion Tip
1 x No.8 Standard Probe/CPITN-C Probe

Periodontal 3 Piece Examination Kit with CPITN-C Probe

SI951ST-1 - 100 kits per box

1 x Front Face Mirror with Percussion Tip
1 x No.8 Standard Probe/CPITN-C Probe
1 x Tweezers



Premium Periodontal Examination Kit with CPITN-C Probe

SI632ST-1 - 100 kits per box

1 x Front Face Mirror with Percussion Tip
1 x No.8 Standard Probe/CPITN-C Probe
1 x London College Tweezers

Complete Periodontal Examination Kit with CPITN-C Probe

SI950ST-1 - 25 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x No.8 Standard Probe/CPITN-C Probe
- 1 x Tweezers
- 2 x Type IIR Face Mask with Ear Loops
- 1 x Patient Bib with Ties
- 2 x Cotton Roll No.2
- 1 x Tissue Wipe
- 1 x Tray (19.5cm x 12.5cm x 2cm)



Suture Removal Kit

SI770ST-1 - 25 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x No.8 Standard Probe/CPITN-C Probe
- 1 x Iris Scissors Curved 11.5cm
- 1 x Iris Toothed Forceps 10cm
- 1 x Tweezers
- 2 x Type IIR Face Mask with Ear Loops
- 2 x Non-Woven Swab 4ply (5cm x 5cm)
- 1 x Tray (19.5cm x 12.5cm x 2cm)



PERIO PROBE/SICKLE EXAMINATION KITS



BEST SELLER

Periodontal 2 Piece Examination Kit

SI552ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle



Periodontal 3 Piece Examination Kit

SI551ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle
- 1 x Tweezers



Premium Periodontal 3 Piece Examination Kit

SI631ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle
- 1 x London College Tweezers

Complete Periodontal Examination Kit

SI550ST-1 - 25 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle
- 1 x Tweezers
- 2 x Type IIR Face Mask with Ear Loops
- 1 x Patient Bib with Ties
- 2 x Cotton Roll No.2
- 1 x Tissue Wipe
- 1 x Tray (19.5cm x 12.5cm x 2cm)



STERILE, SINGLE-USE RESTORATIVE KITS




Labour costs make up a significant portion of a clinic's annual budget with sterilisation of equipment a major contributor to the daily workload. *This is where our kits can make a difference.*

Save time in the steri room

+ Save money on equipment and supplies

+ Improve infection control

Our *restorative instrument* handles are moulded around the stainless steel core, ensuring the instruments display strength when required.



2 Piece Restorative Mixing Kit

SI622ST-1
150 kits per box

- 1 x Spatula DE
- 1 x Dycal Applicator



3 Piece Composite Restorative Kit

SI630ST-1
100 kits per box

- 1 x Flat Plastic DE
- 1 x Ball Burnisher DE
- 1 x Medium Spoon Excavator DE (1.8mm)



5 Piece Amalgam Restorative Kit

SI650ST-1 - 60 kits per box

- 1 x Flat Plastic DE
- 1 x Ball Burnisher DE
- 1 x Medium Spoon Excavator DE (1.8mm)
- 1 x Carver DE
- 1 x Amalgam Plugger Large/Small



6 Piece Composite Restorative Kit

SI660ST-1 - 50 kits per box

- 1 x Flat Plastic DE
- 1 x Ball Burnisher DE
- 1 x Small Spoon Excavator DE (1.2mm)
- 1 x Medium Spoon Excavator DE (1.8mm)
- 1 x PFI49 Burnisher
- 1 x 1/2 Hollenback DE



6 Piece Comprehensive Restorative Kit

SI661ST-1 - 50 kits per box

- 1 x Flat Plastic DE
- 1 x Ball Burnisher DE
- 1 x Medium Spoon Excavator DE (1.8mm)
- 1 x Carver DE
- 1 x Ladmore Plugger DE
- 1 x Amalgam Plugger Large/Small



MDDI - HEAD OFFICE

E: sales@mddiglobal.com | P: +61 (0)7 5453 4433

Unit 1, 24 Hoopers Road, Kunda Park Queensland 4556 AUSTRALIA