



SINGLE-USE INSTRUMENTS CATALOGUE

mddiglobal.com

mddi | 
STERILE | SINGLE USE

CONTENTS

OUR INSTRUMENTS	3
-----------------	---

.....

DIAGNOSTIC INSTRUMENTS	4
------------------------	---

.....

DIAGNOSTIC KITS - CPITN-C PROBE	6
---------------------------------	---

.....

DIAGNOSTIC KITS - PERIO PROBE/SICKLE	7
--------------------------------------	---

.....

RESTORATIVE INSTRUMENTS	9
-------------------------	---

.....

RESTORATIVE KITS	10
------------------	----

.....





Our Instruments are:

- Quick and easy to set up and use
- Prewrapped in a range of convenient kits and instruments
- Sterile at the point of use
- Convenient to dispose of in sharps container
- Reducing the risk of sharps injuries
- Sharp where needed, ensuring the highest quality of treatment for every procedure
- Implemented to significantly reduce the risk of cross-contamination
- GS1 compliant with removable tracking labels
- CE certified and TGA registered:
ARTG 271033, ARTG 272619, ARTG 221241
- 5 year expiry from point of manufacture

MDDI Can Increase Revenue by:

- Minimise decontamination cycles per shift - reduces expenditure on consumables, test kits and utilities
- Putting the RDN back chairside to assist the dentist and patient
- Increasing the number of patients seen in each clinical session
- Improving clinical efficiencies allowing for longer treatment plans and procedures
- Improving the patient journey



STERILE, SINGLE-USE DIAGNOSTIC KITS



Our MDDI diagnostic kits were adopted into the following sectors more than 15 years ago. We have continued to expand our range due to the ongoing popularity of our instruments.

School
Clinics

Government
Hospitals

Flying Doctors
Service

Mobile
Dental Units




Our *diagnostic instruments* are high quality, precision made instruments. The probes and London College Tweezers are manufactured from stainless steel and all of our handles are ABS plastic.

INDIVIDUAL DIAGNOSTIC

SI053ST-1

Front Face Mirror with Percussion Tip

180 per box



SI182ST-1

London College Tweezers

200 per box



SI081ST-1

Tweezers

200 per box





2 Piece Periodontal Examination Kit with CPITN-C Probe

SI952ST-1 - 100 kits per box

1 x Front Face Mirror with Percussion Tip
1 x No. 8 Explorer & CPITN-C Probe

3 Piece Periodontal Examination Kit with CPITN-C Probe

SI951ST-1 - 100 kits per box

1 x Front Face Mirror with Percussion Tip
1 x No. 8 Explorer & CPITN-C Probe
1 x Tweezers



Complete Periodontal Examination Kit with CPITN-C Probe

SI950ST-1 - 25 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x No. 8 Explorer & CPITN-C Probe
- 1 x Tweezers
- 2 x Face Mask with Ear Loops
- 1 x Patient Bib with Ties
- 2 x Cotton Roll No. 2
- 1 x Tissue Wipe
- 1 x Tray (19.5 x 12.5 x 2cm)



Suture Removal Kit

SI770ST-1 - 25 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x No. 8 Explorer & CPITN-C Probe
- 1 x Tweezers
- 2 x Face Mask with Ear Loops
- 2 x Non-Woven Swab 4ply (5 x 5cm)
- 1 x Iris Scissors Curved 11.5cm
- 1 x Iris Toothed Forceps 10cm
- 1 x Tray (19.5 x 12.5 x 2cm)



PERIO PROBE/SICKLE EXAMINATION KITS



BEST SELLER

2 Piece Periodontal Examination Kit

SI552ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Sickle Explorer & Perio Probe



3 Piece Periodontal Examination Kit

SI551ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Sickle Explorer & Perio Probe
- 1 x Tweezers



3 Piece Premium Periodontal Examination Kit

SI631ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Sickle Explorer & Perio Probe
- 1 x London College Tweezers

Complete Periodontal Examination Kit

SI550ST-1 - 25 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Sickle Explorer & Perio Probe
- 1 x Tweezers
- 2 x Face Mask with Ear Loops
- 1 x Patient Bib with Ties
- 2 x Cotton Roll No.2
- 1 x Tissue Wipe
- 1 x Tray (19.5 x 12.5 x 2cm)



Sterile, Single-Use RESTORATIVE Kits




Labour costs make up a significant portion of a clinic's annual budget with sterilisation of equipment a major contributor to the daily workload. *This is where our kits can make a difference.*

Save time in the steri room

Save money on equipment and supplies

Improve infection control

Our restorative instrument handles are moulded around the stainless steel core, ensuring the instruments display strength when required.



3 Piece Restorative Kit

SI630ST-1
100 kits per box

- 1 x Flat Plastic
- 1 x Ball Burnisher
- 1 x Medium Spoon Excavator (1.8mm)



4 Piece Restorative Kit

SI642ST-1 - 80 kits per box


- 1 x Flat Plastic
- 1 x Ball Burnisher
- 1 x 1/2 Hollenback
- 1 x Small Spoon Excavator (1.2mm)



5 Piece Restorative Kit

SI650ST-1 - 60 kits per box

- 1 x Flat Plastic
- 1 x Ball Burnisher
- 1 x Medium Spoon Excavator (1.8mm)
- 1 x Carver
- 1 x Amalgam Plugger Large & Small



6 Piece Restorative Kit

SI661ST-1 - 50 kits per box

- 1 x Flat Plastic
- 1 x Ball Burnisher
- 1 x Medium Spoon Excavator (1.8mm)
- 1 x Carver
- 1 x Ladmore Plugger
- 1 x Amalgam Plugger Large & Small



Our valued distributors



We are constantly looking for new ways to collaborate,
please do not hesitate to get in touch.

MDDI

E: sales@mddiglobal.com | P: 07 5453 4433

Unit 1, 24 Hoopers Road, Kunda Park Queensland 4556 AUSTRALIA



mddiglobal.com